

NANOVAK

NANOVAK® offers global products

Cryo-Pump Large Vacuum Chamber



Hacettepe Üniversitesi Teknokenti, 1596. Cadde No: 95/A, 4. Ar-Ge Binası Ofis 4-5, 06800 Beytepe/Ankara/Turkey
Phone: +90 (312) 299 23 93, Fax: +90 (312) 299 23 94 www.nanovak.com info@nanovak.com



NANOVAK

NANOVAK® offers global products

Cryo-Pump Large Vacuum Chamber

- Cylindrical vacuum chamber made of SS304 material, with electro chemically cleaned and polished surfaces, Φ60 cm x L110 cm in dimensions
- Viton O-rings or copper gaskets on all entries which are durable up to 200 °C
- Φ100 mm watching windows on both door and tank surface
- CF40, CF63, ISO100, ISO250 type feedthroughs on demand
- Wide range (1000 10⁻¹⁰ Torr) vacuum control and measurement system
- Different pumping options, 3500 L/hr cryo-pump, 10⁻⁸ Torr base vacuum level
- Gate valve for isolating cryo pump from tank, automatic gate valve
- PLC controlled vacuum process, automatic (pumping/stopping/venting) control module
- All electronics are seated on sliding 19" drawers,
- 150x110x230cm carrying table on lockable wheels
- 380VAC-16A three phase power, protection against power failures

Hacettepe Üniversitesi Teknokenti, 1596. Cadde No: 95/A, 4. Ar-Ge Binası Ofis 4-5, 06800 Beytepe/Ankara/Turkey Phone: +90 (312) 299 23 93, Fax: +90 (312) 299 23 94 www.nanovak.com info@nanovak.com